



# When you need pure and untouched human CD8<sup>+</sup> T cells

## Dynabeads<sup>®</sup> Untouched<sup>™</sup> Human CD8 T Cells

- Truly untouched cells—no columns required
- High purity, yield, and viability
- Compatible with flow cytometry

The Dynabeads<sup>®</sup> Untouched<sup>™</sup> Human CD8 T Cells kit is the product of choice when you need truly untouched cells for your experiments. This easy-to-use kit is based on negative isolation of CD8<sup>+</sup> T cells by Dynabeads<sup>®</sup> magnetic separation technology. The gentle separation method ensures that your cells are not exposed to the stress of being passed through a column.


### Truly untouched cells

The Dynabeads<sup>®</sup> are used to deplete all the unwanted cells (Figure 1), thus leaving the target cells in the sample untouched and more like the *in vivo* state. The resulting pure, viable, and untouched CD8<sup>+</sup> T cells can be directly analyzed in a flow cytometer and used in any functional assay.


### High purity, yield, and viability

The isolated untouched CD8<sup>+</sup> cells are of very high purity (>85%) and viability (98%) (Figure 2). The cells are readily activated in culture with Dynabeads<sup>®</sup> CD3/CD28 T Cell Expander (Figure 3) and show antigen-specific expansion in culture (Figure 4).

The isolated CD8<sup>+</sup> T cells can be used in any cell-based assays, such as the study of CD8<sup>+</sup> T cell proliferation, apoptosis, and induction of anergy, study of antigen-specific T cells and regulation of CD8<sup>+</sup> T cell cytokine expression, and flow cytometry or fluorescence-activated cell sorting.




**Figure 1**—The gentle and stress-free Dynabeads<sup>®</sup> technology isolates untouched CD8<sup>+</sup> T cells from human peripheral blood mononuclear cells (PBMC). An antibody mix is added to bind to non-CD8<sup>+</sup> T cells (B cells, NK cells, monocytes, platelets, dendritic cells, CD4<sup>+</sup> T cells, granulocytes, and erythrocytes). Dynabeads<sup>®</sup> then bind to these antibody-labeled cells. With the aid of a DynaMag<sup>™</sup> magnet, the bead-bound cells are captured and discarded. The remaining untouched human CD8<sup>+</sup> T cells are ready for flow analysis and any downstream application.







# Negative Cell Isolation



**Figure 2—High cell purity and viability.** CD8<sup>+</sup> cells isolated using the Dynabeads® Untouched™ Human CD8 T Cells kit show high purity (A) and viability (B). Viability was measured using annexin V and propidium iodide staining.



**Figure 3—Activation and proliferation of CD3/CD28–stimulated CD8<sup>+</sup> T cells.** Isolated CD8<sup>+</sup> cells were stimulated in culture with Dynabeads® CD3/CD28 T Cell Expander. Expression of the early activation markers CD25 and CD69 was measured by flow cytometry after 18 hours. After 4 days in culture, 85% of the cells had proliferated, as measured by carboxyfluorescein diacetate succinimidyl ester (CFSE) staining (A). Unstimulated CD8<sup>+</sup> T cells did not proliferate (B).



**Figure 4—Antigen-specific expansion in culture.** Isolated CD8<sup>+</sup> T cells were stimulated with CMV-specific peptide and irradiated autologous peripheral blood mononuclear cells (PBMC). After 11 days, expansion of antigen-specific CD8<sup>+</sup> T cells was analyzed by pentamer staining in a flow cytometer.

## Ordering information

Product	Quantity	Cat. no.
Dynabeads® Untouched™ Human CD8 T Cells*	Processes up to 1 x 10 <sup>9</sup> cells	113-48D
<b>Related products</b>		
Dynabeads® CD3/CD28 T Cell Expander	2 ml	111-31D
CellTrace™ CFSE Cell Proliferation Kit	1 kit	C34554
Vybrant® Apoptosis Assay Kit #2	50 assays	V13241

\* This kit has replaced the Dynabeads® MyPure™ CD8 Kit 2 (Cat. no. 113-41D) and the DYNAL® CD8 Negative Isolation Kit (Cat. no. 113-19D).

For the complete line of fluorescent conjugated antibodies, visit [www.invitrogen.com/antibodies](http://www.invitrogen.com/antibodies). For current prices, visit [www.invitrogen.com](http://www.invitrogen.com). For information on DynaMag™ magnets, visit [www.invitrogen.com/magnets](http://www.invitrogen.com/magnets).

Learn more about the optimal starting point for your T cell research at [www.invitrogen.com/cellisolation](http://www.invitrogen.com/cellisolation).



[www.invitrogen.com](http://www.invitrogen.com)