

PureLink® PCR Purification Kit

**For rapid, efficient purification of PCR
products**

Catalog nos. K3100-01, K3100-02

Rev. date: 20 September 2010
Manual part no. 25-0715

MAN0000446

Contents

| | |
|---|-----------|
| Experienced Users Procedure..... | v |
| Kit Contents and Storage | vi |
| Introduction..... | 1 |
| Product Overview | 1 |
| Experimental Overview | 3 |
| Methods | 5 |
| Purification Procedure..... | 5 |
| Analyzing DNA Yield and Primer Removal..... | 9 |
| Troubleshooting | 11 |
| Appendix | 12 |
| Accessory Products..... | 12 |
| Technical Support | 13 |
| Notes..... | 15 |

Experienced Users Procedure

Note

This quick reference sheet is included for experienced users of the PureLink® PCR Purification Kit. If you are a first time user, follow the detailed protocol provided in this manual.

Purification Procedure

The purification procedure is designed for purifying up to 40 µg dsDNA in a total time of 10–12 minutes using a microcentrifuge capable of centrifuging $>10,000 \times g$.

1. Add 4 volumes of the appropriate PureLink® Binding Buffer with isopropanol (page 7) to 1 volume of PCR (50–100 µL). Mix well.
 2. **Add** sample from Step 1 to a PureLink® Spin Column in a collection tube.
 3. Centrifuge the column at $10,000 \times g$ for 1 minute. dsDNA is bound to the column. Discard the flow through.
 4. **Wash** the column with 650 µL of Wash Buffer with ethanol (page 6).
 5. Centrifuge the column at $10,000 \times g$ for 1 minute. Discard the flow through.
 6. Centrifuge the column at maximum speed for 2–3 minutes to remove any residual wash buffer.
 7. Place the spin column in a clean 1.7-mL elution tube supplied with the kit.
 8. **Elute** with 50 µL of Elution Buffer (add buffer to the center of the column).
 9. Incubate the column at room temperature for 1 minute.
 10. Centrifuge the column at maximum speed for 2 minutes.
 11. The elution tube contains your purified PCR product. Discard the column. Store the purified PCR product at -20°C or use PCR product for the desired downstream application.
-

Kit Contents and Storage

Shipping and Storage All components of the PureLink® PCR Purification Kit are shipped at room temperature. Upon receipt, store all components at room temperature.

Kit Contents The components included in the PureLink® PCR Purification Kit are listed below.
Sufficient reagents are provided in the kit to perform 50 (Cat. no. K3100-01) or 250 (Cat. no. K3100-02) reactions.

| Component | K3100-01 | K3100-02 |
|--|----------|----------|
| Binding Buffer (B2) | 15 mL | 72 mL |
| Binding Buffer High-Cutoff (B3) | 23 mL | 109 mL |
| Wash Buffer (W1) | 16 mL | 80 mL |
| Elution Buffer; 10 mM Tris-HCl, pH 8.5 (E1) | 15 mL | 15 mL |
| PureLink® PCR Spin Columns with Collection Tubes | 50 | 5 × 50 |
| PureLink® Elution Tubes (1.7 mL) | 50 | 5 × 50 |

Intended Use For research use only. Not intended for any animal or human therapeutic or diagnostic use.

Introduction

Product Overview

Description of the System

The PureLink® PCR Purification Kit is designed for rapid and efficient purification of PCR products.

The kit is designed to efficiently remove primers, dNTPs, enzymes, and salts from PCR products in less than 15 minutes. Using the kit with Binding Buffer HC (High-Cutoff) efficiently removes primer dimers or short spurious PCR products (see **PureLink® Binding Buffers**, below).

The purified PCR product is suitable for automated fluorescent DNA sequencing, restriction enzyme digestion, and cloning.

The PureLink® PCR Purification Kit is based on the selective binding of dsDNA to silica-based membrane in the presence of chaotropic salts.

The PCR product is mixed with Binding Buffer to adjust conditions for subsequent dsDNA binding to the PureLink® Spin Column. The dsDNA binds to the silica-based membrane in the column and impurities are removed by thorough washing with Wash Buffer. The dsDNA is then eluted in low salt Elution Buffer or water.

PureLink® Binding Buffers

The PureLink® PCR Purification Kits are supplied with two proprietary buffers:

- Binding Buffer: For routine purifications of 100 bp to 12 kb dsDNA PCR fragments
- Binding Buffer HC (High-Cutoff): For removal of primer dimers or short failed PCR products (<300 bp), eliminating the need for tedious gel purification.

Note: Using Binding Buffer HC reduces the recovery of dsDNA fragments between 300–600 bp, while dsDNA fragments <300 bp do not bind to the PureLink® Spin Column.

Continued on next page

Product Overview, Continued

Advantages


The advantages of using PureLink® PCR Purification Kit are:

- Efficient removal of primers, dNTPs, salts, and enzymes without the need to perform ethanol precipitation
 - Designed to purify PCR products in less than 15 minutes
 - Choice of Binding Buffers for routine purification or selective removal of primer dimers (<300 bp) and short spurious PCR products
 - Reliable performance of the purified PCR products in downstream applications
-

Experimental Overview

Scheme

The flow chart for purifying PCR products using the PureLink® PCR Purification Kit is shown below.



Product Specifications

PureLink[®] PCR Spin Column

| | |
|---------------------------|--|
| Binding Capacity: | 40 µg dsDNA |
| Column Reservoir Capacity | 800 µL |
| Elution Tube Capacity | 1.7 mL |
| Centrifuge Compatibility | Capable of centrifuging >10,000 × <i>g</i> |

System Specifications

| | |
|---------------------------------------|--|
| Starting Material | 50–100 µL PCR product (50 ng–40 µg dsDNA) |
| Elution Volume | 50 µL |
| Separation Range: (Binding Buffer) | 0.1–12 kb from 10–40 mer primers |
| Separation Range: (Binding Buffer HC) | >600 bp from <300 bp PCR fragments and 10–40 mer primers |
| DNA Recovery: | >80% |
| Primer Removal: | >99% |

Methods

Purification Procedure

Introduction The purification procedure is designed for purifying up to **40 µg dsDNA** using a centrifuge in a total time of **10–12 minutes**.

Materials Needed

- 100% isopropanol
 - 96–100% ethanol
 - Binding Buffer (supplied with the kit)
 - Wash Buffer (supplied with the kit)
 - Elution Buffer (supplied with the kit)
 - PureLink® PCR Spin Column and Collection Tubes (supplied with the kit)
 - PureLink® Elution Tubes (supplied with the kit)
 - Sterile, distilled water (pH>7.0)
 - Microcentrifuge capable of achieving >10,000 × g
-



The PureLink® PCR Purification Kit buffers contain guanidine hydrochloride and isopropanol. Always wear a laboratory coat, disposable gloves, and eye protection when handling buffers.

Do not add bleach or acidic solutions directly to solutions containing guanidine hydrochloride or sample preparation waste as it forms reactive compounds and toxic gases when mixed with bleach or acids.

Continued on next page

Purification Procedure, Continued



Follow the recommendations below to obtain the best results:

- Recommended PCR volume is 50–100 μ L
 - Save an aliquot of PCR products before purification to verify and check amplicon on the gel
 - Perform all centrifugation steps at room temperature
 - Pipet the Elution Buffer in the center of the column and perform a 1 minute incubation
 - Always use sterile water with pH 7–8.5, if you are using water for elution
-

Before Starting

- Add 10 mL 100% isopropanol to 15 mL Binding Buffer included with Cat. no. K3100-01.
Add 48 mL 100% isopropanol to 72 mL Binding Buffer included with Cat. no. K3100-02
Store the Binding Buffer with isopropanol at room temperature.
 - Add 2.3 mL isopropanol to 23 mL Binding Buffer HC included with Cat. no. K3100-01.
Add 11 mL 100% isopropanol to 109 mL Binding Buffer HC included with Cat. no. K3100-02
Store the Binding Buffer HC with isopropanol at room temperature.
 - Add 64 mL 96–100% ethanol to 16 mL Wash Buffer included with Cat. no. K3100-01.
Add 320 mL 96–100% ethanol to 80 mL Wash Buffer included with Cat. no. K3100-02
Store the Wash Buffer with ethanol at room temperature.
-

Continued on next page

Purification Procedure, Continued

Binding DNA

1. Add 4 volumes of PureLink® Binding Buffer with isopropanol (see page 6) or Binding Buffer HC with isopropanol (see page 6) to 1 volume of the PCR product (50–100 µL). Mix well.
 2. Remove a PureLink® Spin Column in a Collection Tube from the package.
 3. Add sample with appropriate Binding Buffer from Step 1 to the PureLink® Spin Column.
 4. Centrifuge the column at room temperature at $10,000 \times g$ for 1 minute.
 5. Discard the flow through and place the spin column into the collection tube.
 6. Proceed to **Washing DNA**, below.
-

Washing DNA

1. Add 650 µL of Wash Buffer with ethanol (page 6) to the column.
 2. Centrifuge the column at room temperature at $10,000 \times g$ for 1 minute. Discard the flow through from the collection tube and place the column into the tube.
 3. Centrifuge the column at maximum speed at room temperature for 2–3 minutes to remove any residual Wash Buffer. Discard the collection tube.
 4. Proceed to Eluting DNA, next page.
-

Continued on next page

Purification Procedure, Continued

Eluting DNA

1. Place the spin column in a clean 1.7-mL PureLink® Elution Tube supplied with the kit.
2. Add 50 μ L of Elution Buffer (10 mM Tris-HCl, pH 8.5) or sterile, distilled water (pH >7.0) to the center of the column.
3. Incubate the column at room temperature for 1 minute.
4. Centrifuge the column at maximum speed for 2 minutes.
5. The elution tube contains your purified PCR product. Remove and discard the column. The recovered elution volume is ~48 μ L.
6. Store the purified PCR product at -20°C or use PCR product for the desired downstream application.

Examples of efficient primer removal using the different Binding Buffers are described on page 10.

Analyzing DNA Yield and Primer Removal

DNA Yield

After purification with PureLink® PCR Purification Kit, the yield of purified dsDNA can be estimated by agarose gel electrophoresis or Quant-iT™ DNA Assay Kits.

Agarose Gel Electrophoresis

To estimate the yield, perform agarose gel electrophoresis of the purified PCR product and known quantities of DNA fragment of the same size. Compare the band intensity of the purified PCR product with the standard DNA fragments.

Quant-iT™ DNA Assay Kits

The Quant-iT™ DNA Assay Kits (see page 12 for ordering information) provide a rapid, sensitive, and specific method for dsDNA quantitation with minimal interference from RNA, protein, ssDNA (primers), or other common contaminants that affect UV absorbance.

The kit contains a state-of-the-art quantitation reagent, pre-diluted standards for standard curve, and a ready-to-use buffer. The assay is performed in a microtiter plate format and is designed for reading in standard fluorescent microplate readers. Follow manufacturer's recommendations to perform the assay.

Primer Removal

The efficiency of primer removal can be estimated by agarose gel electrophoresis as described in the examples shown on the next page.

The WAVE® System is an ideal method to estimate the efficiency of primer removal. The WAVE® System is an automated DHPLC (denatured high-performance liquid chromatography) system.


Continued on next page

Analyzing DNA Yield and Primer Removal, Continued

Example with Binding Buffer

An example of efficient primer removal using Binding Buffer is shown below.


A mixture of 100 bp DNA Ladder (see page 12) with an excess of a 37 mer primer was purified using Binding Buffer with the PureLink[®] PCR Purification Kit as described in the manual. The mixture was analyzed using agarose gel electrophoresis before purification (lane 1) and after purification (lane 2).



Example with Binding Buffer HC

An example of efficient removal of <300 bp fragments using Binding Buffer HC is shown below.

PCR product was purified using Binding Buffer HC with PureLink[®] PCR Purification Kit as described in the manual. The PCR product was analyzed by agarose gel electrophoresis and shows efficient purification of a 965-bp PCR product from a 232-bp PCR fragment and primers. The figure below shows a DNA Ladder (lane 1), the PCR product before purification (lane 2) and the PCR product after purification (lane 3).



Troubleshooting

| Observation | Cause | Solution |
|--|-------------------------------------|--|
| Low DNA yield | PCR conditions not optimized | Check amplicon on gel to verify the PCR product prior to purification. |
| | Incorrect binding conditions | For efficient DNA binding always mix 1 volume of PCR (50–100 μ L) with 4 volumes of Binding Buffer. Be sure to add 100% isopropanol to the Binding Buffer as described on page 7. |
| | Ethanol not added to Wash Buffer | Be sure to add 96–100% ethanol to Wash Buffer as described on page 6. |
| | Incorrect elution conditions | Add elution buffer to the center of the column and perform incubation for 1 minute with elution buffer before centrifugation. |
| Presence of primer dimers | Incorrect Binding Buffer used | For efficient removal of primer dimers or short failed PCR products (<300 bp), use Binding Buffer HC. This buffer is specifically designed to remove <300 bp DNA fragments eliminating the need for gel purification. |
| Inhibition of downstream enzymatic reactions | Presence of ethanol in purified DNA | Traces of ethanol from the Wash Buffer can inhibit downstream enzymatic reactions. To remove Wash Buffer, discard Wash Buffer flow through from the collection tube. Place the spin column into the collection tube and centrifuge the spin column at maximum speed for 2–3 minutes to completely dry the column. |

Appendix

Accessory Products

Introduction The following products may be used with the PureLink® PCR Purification Kit. For details, visit www.invitrogen.com or contact **Technical Support** (see page 13).

| Item | Amount | Cat. no. |
|---|------------------|-----------|
| PureLink® 96 PCR Purification Kit | 4 × 96 reactions | K3100-96 |
| Platinum® <i>Taq</i> DNA Polymerase High Fidelity | 100 reactions | 11304-011 |
| Platinum® <i>Taq</i> DNA Polymerase | 100 reactions | 10966-018 |
| UltraPure™ DNase/RNase-free Distilled Water | 500 mL | 10977-015 |
| Quant-iT™ DNA Assay Kit, High Sensitivity | 1000 assays | Q33120 |
| Quant-iT™ DNA Assay Kit, Broad-Range | 1000 assays | Q33130 |
| PureLink™ 96 Receiver Plate | 50 | 12193-025 |
| 100 bp DNA Ladder | 50 µg | 15628-019 |

Technical Support

Web Resources



Visit the Invitrogen website at www.invitrogen.com for:

- Technical resources, including manuals, vector maps and sequences, application notes, SDSs, FAQs, formulations, citations, handbooks, etc.
- Complete technical support contact information
- Access to the Invitrogen Online Catalog
- Additional product information and special offers

| | | |
|---|---|--|
| Corporate Headquarters: 5791 Van Allen Way Carlsbad, CA 92008 USA Tel: 1 760 603 7200 Tel (Toll Free): 1 800 955 6288 Fax: 1 760 602 6500 E-mail: tech_support@invitrogen.com | Japanese Headquarters: LOOP-X Bldg. 6F 3-9-15, Kaigan Minato-ku, Tokyo 108-0022 Tel: 81 3 5730 6509 Fax: 81 3 5730 6519 E-mail: jpinfo@invitrogen.com | European Headquarters: Inchinnan Business Park 3 Fountain Drive Paisley PA4 9RF, UK Tel: +44 (0) 141 814 6100 Tech Fax: +44 (0) 141 814 6117 E-mail: eurotech@invitrogen.com |
|---|---|--|

SDS Information

SDSs (Safety Data Sheets) are available on our website at www.invitrogen.com/sds.

Continued on next page

Technical Support, Continued

Limited Warranty

Invitrogen (a part of Life Technologies Corporation) is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Support Representatives.

All Invitrogen products are warranted to perform according to specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits the Company's liability to only the price of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. The Company reserves the right to select the method(s) used to analyze a product unless the Company agrees to a specified method in writing prior to acceptance of the order.

The Company makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore the Company makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Support Representatives.

Life Technologies Corporation shall have no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

©2010 Life Technologies Corporation. All rights reserved.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

Notes

Notes



Corporate Headquarters

5791 Van Allen Way
Carlsbad, CA 92008

T: 1 760 603 7200

F: 1 760 602 6500

E: tech_support@invitrogen.com

For country-specific contact information visit our web site at www.invitrogen.com