

# TOPO® TA Cloning® Kits

Publication Part No. 25-0185

MAN0005242

Rev. Date: 13 March 2012

## Produce PCR products

1. Produce PCR products using *Taq* polymerase. Include a final extension step of 7–10 minutes to ensure 3' adenylation.
2. Verify there is only one PCR product of the proper size in each reaction by agarose gel electrophoresis.

## Perform TOPO® cloning

1. Set up a 6  $\mu$ L TOPO® Cloning reaction for each PCR product:

Reagent	Volume
Fresh PCR product	0.5–4 $\mu$ L
Salt solution	1 $\mu$ L
Water	add to a final volume of 5 $\mu$ L
TOPO® Vector	1 $\mu$ L

2. Mix gently and incubate for 5 minutes at room temperature.
3. Place tubes on ice. Proceed to **Transform One Shot® cells** or **Electroporate One Shot® cells**.

## Transform One Shot® cells

1. Thaw One Shot® *E. coli* or equivalent competent cells on ice.
2. Add 2  $\mu$ L of the TOPO® Cloning reaction to a vial of One Shot® *E. coli* or equivalent competent cells and mix gently.
3. Incubate on ice for 5–30 minutes.
4. Heat-shock the cells for 30 seconds at 42°C without shaking. Proceed to **Recover and plate cells**.

## Electroporate One Shot® cells

1. Thaw competent cells on ice.
  2. Dilute the TOPO® Cloning reaction 4-fold with sterile water.
  3. Add 2  $\mu\text{L}$  of the diluted TOPO® Cloning reaction to a vial of One Shot® Electrocomp™ *E. coli* or equivalent competent cells and mix gently.
  4. Electroporate your samples. Refer to the manual if arcing occurs.
  5. Proceed to **Recover and plate cells**.
- 

## Recover and plate cells

1. Add 250  $\mu\text{L}$  of room temperature S.O.C. Medium to the cells. For electroporated cells, transfer the sample to a 15-mL culture tube.
  2. Cap the tubes and shake at 37°C for 1 hour.
  3. Spread 10–50  $\mu\text{L}$  from each transformation on prewarmed LB plates containing X-gal, and either 50–100  $\mu\text{g}/\text{mL}$  ampicillin or 50  $\mu\text{g}/\text{mL}$  kanamycin. For TOP10F' *E. coli*, add both IPTG and X-gal to the LB plates.
  4. Incubate plates overnight at 37°C.
  5. Pick ~10 white or light blue colonies for analysis. Do not pick dark blue colonies. Proceed to **Analyze positive clones**.
- 

## Analyze positive clones

1. Culture 10 colonies in LB medium containing 50–100  $\mu\text{g}/\text{mL}$  ampicillin or 50  $\mu\text{g}/\text{mL}$  kanamycin.
  2. Isolate plasmid DNA using your method of choice.
  3. Analyze the plasmid by restriction analysis or by sequencing to confirm the presence and correct orientation of the insert.
- 

For more information, refer to [www.lifetechnologies.com/manuals](http://www.lifetechnologies.com/manuals)

© 2012 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

For support visit [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support)  
[www.lifetechnologies.com](http://www.lifetechnologies.com)

