



Comments for pBudCE4.1
4595 nucleotides

CMV promoter: bases 7-594

CMV Forward priming site: bases 544-564

T7 promoter/priming site: bases 638-657

CMV multiple cloning site: bases 664-713

myc epitope: bases 719-748

6xHis tag: bases 764-782

SV40 polyadenylation sequence: bases 803-933

Zeocin™ resistance gene: bases 1063-1437 (complementary strand)

EM7 promoter: bases 1456-1510 (complementary strand)

SV40 early promoter: bases 1547-1869 (complementary strand)

EF-1α promoter: bases 1885-3051

EF-1α Forward priming site: bases 2999-3019

EF-1α multiple cloning site: bases 3062-3126

V5 epitope: bases 3127-3168

6xHis tag: bases 3178-3195

BGH Reverse priming site: bases 3218-3235 (complementary strand)

BGH polyadenylation sequence: bases 3224-3447

pUC origin: bases 3521-4194