



**Comments for:
(no insert)**

T7 promoter/priming site:
 Sp6 promoter/priming site (c):
 CMV promoter (c):
lacUV5 promoter (c):
 pUC origin:
 Ampicillin (*bla*) resistance gene (c):
bla promoter (c):
 f1 origin (c):
 SV40 early polyadenylation signal (c):
 (c) = complementary strand

pCMV•SPORT1
 4380 nucleotides

16-35
 156-175
 262-855
 856-883
 1123-1821
 2146-3006
 3001-3058
 3294-3749
 3876-4142

pCMV•SPORT2
 4475 nucleotides

11-30
 151-170
 252-845
 900-927
 1335-2004
 2191-3051
 3046-3103
 3339-3794
 3921-4187

pCMV•SPORT4
 4387 nucleotides

11-30
 151-169
 252-845
 --
 1243-1916
 2103-2963
 2959-3015
 3251-3706
 3833-4099