



**Comments for pIND**  
**5023 nucleotides**

Ecdysone/glucocorticoid

Response Elements (5 x E/GREs): bases 1-174

Minimal Heat Shock Promoter: bases 181-475

Multiple Cloning Site: bases 476-601

BGH Reverse priming site: bases 603-620

BGH Polyadenylation sequence: bases 602-816

f1 origin: bases 879-1293

SV40 promoter/origin: bases 1357-1681

Neomycin ORF: bases 1717-2511

SV40 Polyadenylation sequence: bases 2527-2766

pUC origin: bases 3198-3871

Ampicillin ORF: bases 4016-4876 (complementary strand)