



**Comments for pLenti6.2/C-Lumio™/V5-GW/lacZ
10363 nucleotides**

RSV enhancer/promoter: bases 1-229
HIV-1 5' LTR: bases 230-410
5' splice donor: base 520
HIV-1 psi (ψ) packaging signal: bases 521-565
HIV-1 Rev response element (RRE): bases 1075-1308
3' splice acceptor: base 1656
3' splice acceptor: base 1684
CMV promoter: bases 1873-2462
attB1 site: bases 2506-2530
lacZ ORF: bases 2582-5633
attB2 site: bases 5648-5672
V5 epitope: bases 5686-5727
Lumio™ tag: bases 5746-5763
PGK promoter: bases 5815-6321
EM7 promoter: bases 6332-6398
Blasticidin resistance gene: bases 6405-6797
ΔU3/HIV-1 3' LTR: bases 6883-7117
 ΔU3: bases 6883-6936
 Truncated HIV-1 3' LTR: bases 6884-7117
SV40 polyadenylation signal: bases 7189-7391
bla promoter: bases 8179-8283
Ampicillin (*bla*) resistance gene: bases 8278-9138
pUC origin: bases 9283-9956