



Comments for pLenti6.3/V5-TOPO
7692 nucleotides

- RSV/5' LTR hybrid promoter: bases 1-410
- RSV promoter: bases 1-229
- HIV-1 5' LTR: bases 230-410
- HIV-1 psi (ψ) packaging signal: bases 521-565
- HIV-1 Rev response element (RRE): bases 1075-1308
- cPPT: bases 1800-1922
- CMV promoter: bases 1934-2518
- TOPO[®] cloning site: bases 2557-2566
- V5 epitope: bases 2629-2670
- WPRE: bases 2689-3286
- SV40 promoter: bases 3297-3605
- EM7 promoter: bases 3660-3726
- Blasticidin resistance gene: bases 3727-4125
- Δ U3/3' LTR: bases 4211-4445
 - Δ U3: bases 4211-4264
 - 3' LTR: bases 4265-4445
- SV40 polyadenylation signal: bases 4517-4648
- bla* promoter: bases 5507-5605
- Ampicillin (*bla*) resistance gene: bases 5606-6466
- pUC origin: bases 6611-7284