



**Comments for pLenti6/TR
8355 nucleotides**

RSV/5' LTR hybrid promoter: bases 1-410

RSV promoter: bases 1-229

HIV-1 5' LTR: bases 230-410

5' splice donor: base 520

HIV-1 psi (ψ) packaging signal: bases 521-565

HIV-1 Rev response element (RRE): bases 1075-1308

3' splice acceptor: base 1656

3' splice acceptor: base 1684

CMV promoter: bases 1809-2392

Rabbit β-globin intron II (IVS): bases 2552-3124

tetR gene: bases 3208-3855

SV40 early promoter and origin: bases 3960-4269

EM7 promoter: bases 4324-4390

Blasticidin resistance gene: bases 4391-4789

ΔU3/3' LTR: bases 4875-5109

ΔU3: bases 4875-4928

3' LTR: bases 4929-5109

SV40 polyadenylation signal: bases 5181-5315

bla promoter: bases 6171-6269

Ampicillin (*bla*) resistance gene: bases 6270-7130

pUC origin: bases 7275-7948